

DNA extractions using the *E.Z.N.A.* kit

Lucas A. Nell

May 26, 2016

Modified instructions specific for extracting DNA from mouse liver samples

Part 1

Before beginning:

- Set heat block to 55°C.
- Fill an ice pan with ice.

Instructions:

1. Transfer ~20 mg tissue to a 1.5 mL microcentrifuge tube. Keep sample in ice pan when not in use.
2. Add 200 μ L TL Buffer.
3. Grind sample using disposable pestle and motor.
4. Add 25 μ L OB Protease Solution. Vortex to mix thoroughly.
5. Incubate at 55°C in the heat block, and vortex the sample every 20-30 minutes (when possible). Allow sample to incubate overnight, or at least 3 hours if you're in a rush.

Part 2

Before beginning (~1 hr before starting):

- Set water bath to 70°C, fill with distilled water, and allow it to reach 70°C.
- Vortex Elution Buffer and transfer enough of it for all samples ($[\# \text{ samples} + 1] \times 100\mu\text{L}$) to microcentrifuge tube(s). Place in heated water bath.
- Dilute HBC Buffer with isopropanol (see bottle for instructions).
- Dilute DNA Wash Buffer with ethanol (see bottle for instructions).

Instructions:

Note: "HiBind DNA Mini Column" and "Collection Tube" are those provided in the *E.Z.N.A.* kit.

1. Remove sample from heat block, vortex to mix, and briefly spin down to remove moisture from inside of the lid.
2. Add 5 μ L RNase A, and let sit at room temperature for 2 minutes.
3. Centrifuge at maximum speed ($\geq 10,000 \times g$) for 5 minutes at room temperature.
4. Transfer the supernatant to a sterile 1.5 mL microcentrifuge tube. Do not disturb or transfer any of the insoluble pellet.

5. Add 220 μ L BL Buffer. Vortex to mix thoroughly. A wispy precipitate may form upon the addition of BL Buffer. This does not interfere with DNA recovery.
6. Incubate at 70°C in water bath for 10 minutes.
7. Add 220 μ L 100% ethanol. Vortex to mix thoroughly.
8. Insert a HiBind DNA Mini Column into a 2 mL Collection Tube.
9. Transfer the entire sample in the microcentrifuge tube to the HiBind DNA Mini Column, including any precipitates that may have formed.
10. Centrifuge at maximum speed for 1 minute. Discard the filtrate; reuse the collection tube.
11. Add 500 μ L HBC Buffer.
12. Centrifuge at maximum speed for 30 seconds. Discard the filtrate *and* collection tube.
13. Insert the HiBind DNA Mini Column into a new 2 mL Collection Tube.
14. DNA wash:
 - a. Add 700 μ L DNA Wash Buffer.
 - b. Centrifuge at maximum speed for 30 seconds.
 - c. Discard the filtrate; reuse the collection tube.
15. *Run a second DNA wash.*
16. Centrifuge the empty HiBind DNA Mini Column at maximum speed for 2 minutes to dry the column. This step is critical for removal of trace ethanol that may interfere with downstream applications.
17. Transfer the HiBind DNA Mini Column into a nuclease-free 1.5 mL microcentrifuge tube.
18. DNA elution:
 - a. Add 50 μ L Elution Buffer heated to 70°C.
 - b. Let sit at room temperature for 2 minutes.
 - c. Centrifuge at maximum speed for 1 minute.
19. *Run a second DNA elution.*
20. Store eluted DNA at -20°C.